

INTRODUCTION

Laboratory animals play a vital role in biomedical research. Researchers use it as a model for developing a cure for different diseases. However, persons involved in the care and use of animals may be exposed to different hazards—allergens, zoonoses, and physical hazards.

Occupational hazard is inevitable during the course of animal research. Without the proper equipment for protection, these hazards may cause serious damage to the personnel in the long run. Esco offers **VIVA® Animal Research Workstations**, providing both personnel and environment protection from exposure to biological hazards and unpleasant odors.

VIVA® Universal Animal Containment Workstation

A containment equipment that provides Biosafety Cabinet Class II.

VIVA® Dual Access Animal Containment Workstation

A workstation with accessible front and back sides that is suitable for cage changing and animal transfer procedure.

VIVA® Bedding Disposal Animal Containment Workstation

A workstation for cage changing and bedding disposal procedures.


List of Applications

- Animal Behavioral Studies
- Animal Breeding
- Drug Development
- Cosmetic Testing
- Biotechnology
- Academic Research



ESCO GLOBAL NETWORK

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD




ESCO
SCIENTIFIC

Esco Animal Research Workstations



Join us on Social Media and Download our Apps!



Filtra Check App
The World's First Chemical
Evaluation App



Esco Social Media
and Apps

Esco Animal Research Workstations



OVERVIEW OF MODELS

VIVA® Universal Animal Containment Workstation

- Sentinel™ Gold Microprocessor Control System
- Ergonomic, ADA-compliant
- Sloped Front Angle
- Available Sizes: 4 and 6 ft



Model: VA2-4A_

VIVA® Dual Access Animal Containment Workstation

- Sentinel™ Gold Microprocessor Control System
- Ergonomic, ADA-compliant
- Advanced Work Tray Design
- Available Sizes: 4 and 5 ft



Model: VDA-4A_

VIVA® Bedding Disposal Animal Containment Workstation

- Sentinel™ Silver Microprocessor Control System
- Integrated Waste Bin
- Nanocarb Activated Carbon Filter for Removing Odor
- Available Size: 4 ft only



Model: VBD-4A_

OPTIONS AND ACCESSORIES



Service Fixtures



Electrical Outlet



Stainless Steel IV Bar
(for VA2 only)



UV Lamp
(for VA2 only)



Side Shield
(for VDA only)



Feed Hopper
(for VDA only)



Foldable Side Tray
(for VDA only)

KEY FEATURES:

- Energy-efficient DC ECM Motor
- ULPA Filter (ISO Class 3 Work Zone)
- Low Noise Level
- ELISA-verified Containment
- Isocide™ Antimicrobial Powder Coating
- Airflow Monitoring System