



A-Series Mechanical Pipette

Transfer with ease and accuracy



A-Series Mechanical Pipette



Ushering in a new age of traceability research, the A-Series Mechanical Pipette is a new line of pipettes aimed at improving ergonomics and comfort while maintaining its accuracy and reliability.

KEY FEATURES

Control Button

- Very low operating force
- Colour-coded based on volume
- Positioned for perfect ergonomics

Volume Adjustment

- Turn to adjust volume

Ejector

- Low tip ejection force of less than or equal to 4 newtons
- Positioned for perfect ergonomics

Volume Display


- 4-digit display
- Magnified for comfortable viewing

Quick Connection Tip

- Remove lower part quickly

Tip Cone

- Spring-loaded tip cone for ease of connecting and ejecting pipette tips



PRODUCT FEATURES

- 4-digit volume display
- Magnified display for comfortable viewing
- Chemical resistant
- Up to 121°C sterilization
- Very low operating force
- Colour-coded according to volume
- Positioned for perfect ergonomics
- Turn to adjust volume
- Antimicrobial surface
- No need for greasing of the pipette piston after autoclaving




Ordering Information

Pipette

ITEM CODE	Model Code	Description	Accuracy	Repeatability
1150488	PIPETTE A2, L-2107	Volume Range of 0.1 - 2 μ l, 1 pcs	\pm 12.0-1.5%	\leq 6.00-0.70%
5171116		Volume Range of 0.1 - 2 μ l, 1 Pack (20 pcs per pack)		
1150489	PIPETTE A10, L-2108	Volume Range of 0.5 - 10 μ l, 1 pcs	\pm 2.5-1.0%	\leq 1.50-0.40%
5171117		Volume Range of 0.5 - 10 μ l, 1 Pack (20 pcs per pack)		
1150490	PIPETTE A20, L-2109	Volume Range of 2 - 20 μ l, 1 pcs	\pm 2.5-1.0%	\leq 1.50-0.30%
5171118		Volume Range of 2 - 20 μ l, 1 Pack (20 pcs per pack)		
1150491	PIPETTE A100, L-2111	Volume Range of 10 - 100 μ l, 1 pcs	\pm 1.8-0.8%	\leq 0.50-0.15%
5171119		Volume Range of 10 - 100 μ l, 1 Pack (20 pcs per pack)		
1150492	PIPETTE A200, L-2112	Volume Range of 20 - 200 μ l, 1 pcs	\pm 1.8-0.6%	\leq 0.50-0.15%
5171120		Volume Range of 20 - 200 μ l, 1 Pack (20 pcs per pack)		
1150493	PIPETTE A1000, L-2113	Volume Range of 100 - 1000 μ l, 1 pcs	\pm 1.5-0.6%	\leq 0.30-0.15%
5171121		Volume Range of 100 - 1000 μ l, 1 Pack (20 pcs per pack)		
1150494	PIPETTE A5000, L-2114	Volume Range of 500 - 5000 μ l, 1 pcs	\pm 1.2-0.6%	\leq 0.30-0.15%
5171122		Volume Range of 500 - 5000 μ l, 1 Pack (20 pcs per pack)		
1150495	PIPETTE A10ML, L-2115	Volume Range of 1 - 10 mL, 1 pcs	\pm 2.5-0.6%	\leq 0.30-0.15%
5171123		Volume Range of 1 - 10 mL, 1 Pack (20 pcs per pack)		
Complete Set				
5171125		Complete Set of 8 Pipettes (A2, A10, A20, A100, A200, A1000, A5000, A10ML) & 1 Pipette Stand (HW-9)		
Accessories				
1150496	PIPETTE STAND HW-9, L-2039A	Pipette stand that holds up to 9 Pipettes		
5171124		Pipette stand that holds up to 9 Pipettes, 1 Pack (20 pcs per pack)		

ESCO LIFESCIENCES GROUP
42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



*Follow us on social media, download our apps,
and scan the QR code for more info.*



@Escolifesciences



@Escolifesciences



@Escolifesci



@Esco



@Escolifesciences



@Escolifesciences



Esco Lifesciences



Esco Lifesciences

ESCO[®]
LIFESCIENCES GROUP

Esco Micro Pte. Ltd. • 19 Changi South Street 1 • Singapore 486 779
Tel +65 6542 0833 • mail@escolifesciences.com
www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA
Tel: +1 215-441-9661 • Fax 484-698-7757
eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania
| Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam