



# Pipette Tips Compatibility Chart

Compatible Pipette Tips			Esco Electronic Pipette										
			Single Channel						Multi-Channel				
			1710001	1710002	1710003	1710004	1710005	1710006	1710007	1710010 & 1710013	1710011 & 1710014	1710012 & 1710015	
Brand	Model	Volume	0.2-10 µL	0.5-20 µL	2-100 µL	5-200 µL	10-300 µL	10-500 µL	50-1000 µL	0.2-10 µL	5-100 µL	10-300 µL	
Oxygen	T-300-RS	10 µL	✓	✓							✓		
	TXL-10	10 µL	✓	✓									
	T-300	10 µL	✓	✓									
	T400	20 µL	✓	✓									
	TF-400	20 µL	✓	✓									
	T-200-Y-R-S	200 µL			✓	✓	✓					✓	✓
	T-200-C	200 µL			✓								
	TF-100	100 µL			✓								
	TF-200	200 µL			✓	✓							
	T-350-C	300 µL						✓					
	TF-350	300 µL						✓					
	T-1000-B-R-S	1000 µL							✓	✓			
T-1000-C	1000 µL							✓	✓				
TF-1000	1000 µL							✓	✓				
QSP	T104RS-Q	10 µL	✓	✓							✓		
	TF140-200-Q	200 µL		✓	✓	✓	✓					✓	✓
	TF112-1000Q	1000 µL							✓	✓			
Eppendorf	0030076141	10 µL	✓	✓							✓		
	0030076150	200 µL		✓	✓	✓	✓					✓	✓
	0030076176	1000 µL							✓	✓			
Sartorius	790011	10 µL	✓	✓							✓		
	790200	200 µL		✓	✓	✓	✓					✓	✓
	791000	1000 µL							✓	✓			
Kirgen	KG1031	10 µL	✓	✓							✓		
	KG5232	200 µL		✓	✓	✓	✓					✓	✓
	KG5333	1000 µL							✓	✓			
Biologix	21-0010	10 µL	✓	✓							✓		
	21-0200	200 µL		✓	✓	✓	✓					✓	✓
	21-1250	1000 µL							✓	✓			

The above data were obtained from in-house tests conducted using tips in forward-mode which conforms to ISO 8655 international standards.

