



The All-in-One Solution For Every Cell Processing Need

A **Cell Processing Isolator (CPI)** is an integrated system that combines several types of equipment into one isolated solution. General bioprocessing equipment such as CO₂ incubator and centrifuge can be integrated in the isolator. These equipment are usually used in a virus/vaccine production facility or cell therapy applications for advanced therapy medicinal products (ATMP).

As a cGMP isolator and an aseptic containment system, it is pressure tested and is capable of biodecontamination procedure. The CPI is designed to isolate the process to ensure operator safety, without compromising the product quality. It also provides the sterile environment (ISO Class 5/Grade A) that is required in carrying out aseptic processes.

It is easily customizable to meet the requirements of each client. Conversion of CPI into an isolator incubator is also available with temperature, relative humidity, CO₂, N₂, and suppressed O₂ control.

Other options include temperature and CO₂ control, Esco cell therapy equipment, or as a stand-alone isolator with third-party equipment.

For more information, please visit www.escopharma.com

Esco GB Ltd – UK

Unit 2 R-Evolution @ Gateway 36, Kestrel Way, Barnsley, S70 5SZ

Tel: +44 (0) 1226 360799 • Email: egb.info@escoglobal.com • Website: www.escopharma.com

Esco Global Offices: Bangladesh | Cameroon | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

ESCO

WORLD CLASS. WORLDWIDE.